
de - Gebrauchsanweisung	5
en - Instructions for use	7
fr - Mode d'emploi	9
it - Istruzioni per l'uso	11
es - Instrucciones de uso	14
pt - Instruções de utilização.....	16
nl - Gebruiksaanwijzing	18
da - Brugsanvisning	21
sv - Bruksanvisning	23
fi - Käyttöohje	25
cs - Návod k použití	27
sk - Návod na použitie	29
sr - Uputstvo za upotrebu	31
hu - Használati útmutató	33
ro - Instrucțiuni de utilizare	36
bg - Инструкции за употреба	38
el - Οδηγίες χρήσης	40

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ENGLISH

Intended purpose

The coaxial cannula is intended for introducing and guiding of compatible needles or cannulas. In this combination it is intended for soft tissue biopsies and placement of marking clips. It is intended for single-use and for transient invasive application.

• Indications

Indications are anomalies in soft tissues as breast, prostate, liver, kidney, lung or thyroid.

• Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder and anticoagulant treatment. In case of increased risks, application depends on decision of healthcare professionals.

• Residual risks and side effects

Residual risks and undesirable side effects are general puncture-related complications such as bleeding, infection, injury to surrounding tissue/organ and pain, as well as pneumothorax and air embolism if used on lung.

• Intended users

Intended users are health care professionals only (physicians).

• Intended patients

Intended patients are adults.

Clinical benefits

Visibility to ultrasound, X-ray (mammography, computer tomography CT) or magnetic resonance imaging MRI.

Allows for multiple tissue removal with only one penetration point and minimal trauma.

Ensures minimal risk of cell displacement for biopsy procedures.

Easy penetration of dense tissue due to the special edge (grinding at tip of coaxial cannula).

Description

The BIP Coaxial Cannula is made of medical-grade stainless steel or titanium (MR-compatible variants) and consists of two parts, an outer cannula with a plastic handle and an inner stylet with a central trocar tip (see Fig. 3). The plastic cap of the stylet is color coded (e.g. CC14[yyy] with 14G stylet – green hub), matching the color code and size of the compatible biopsy needle. Please note that the diameter of the outer cannula is larger than indicated by the color.

The outer cannula is equipped with centimeter depth marks for direct visual assessment of insertion depth, and is equipped with a luer lock connection that enables to connect syringes or tubing to the proximal end of the cannula for optional aspiration (see Fig. 3).

The BIP Coaxial Cannula is sterile, single-packaged in peel pouches and is available in various diameters and lengths (see Fig. 2 and Tab. 1). Sales units are additionally packaged in protective cardboard box.

Compatibility

The BIP Coaxial Cannula is specially designed for use in combination with the BIP EvoCore®, the BIP HistoCore® and the BARD® Magnum® biopsy devices and respective biopsy needles. In addition, BIP marking clips can be applied through the coaxial cannula (see Tab. 1).

The BIP Coaxial Cannula can be used with conventional imaging using ultrasound or X-ray (mammography, computer tomography CT) or magnetic resonance imaging MRI. If the BIP Coaxial Cannula shall be used in magnetic resonance devices (MR systems), ensure that a MR compatible variant is used (see Tab. 1, CC[xx][yyy]MR).

Application

1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. The application normally takes place through healthy, uninjured skin. When using coaxial cannulas with a large diameter, an incision should be made beforehand with a scalpel according to the discretion of the user. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

Prior to use, the gauge size (external diameter) and the length of the BIP Coaxial Cannula should be determined that is required for the particular procedure to be performed (see Tab. 1).

2. Remove the coaxial cannula from sterile package

Check sterile packaging and the product for expiration date and visual damages before use.

Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly  and remove the product aseptically.

3. Positioning

Insert the BIP Coaxial Cannula under the control of suitable compatible imaging procedure until the stylet tip is directly in front of the area to be treated. Hold the plastic handle of the coaxial cannula and turn the inner stylet on the handle by 90° counter-clockwise and pull out (see Fig. 3). The outer cannula remains in the patient and serves as a guide for the biopsy needle or other application cannula.

4. Performing biopsies or placement of marking clips

Insert compatible biopsy needle or application cannula into the guide cannula until stop.

Observe the instructions for use of the respective biopsy device, biopsy needle or marking clip being used.

Perform the biopsy or placement of marking clip and remove the biopsy needle or application cannula from the guide cannula again.

NOTE: When used in MR environment, position the patient outside the 5 Gauss line before inserting the biopsy needle, as biopsy needle is not MR compatible.

When taking multiple samples, leave the guide cannula in position, reinsert the biopsy needle into the guide cannula and perform the next biopsy.

5. Removal and disposal

After procedure remove coaxial cannula from patient and dispose of in accordance with local provisions and safety regulations.

REACH declaration at <https://www.bipmedical.com/en/about-us/regulations/> must be observed.

⚠ Warnings

- For those organs as listed in the indications, safety and performance is demonstrated on sufficient clinical data. If the device is to be used on other similar soft tissues, except those body regions listed as contraindications, safe application depends on decision of healthcare professionals.
- The coaxial cannula must not be used and must be disposed of if the product is damaged or if the expiration date has elapsed or if the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to insertion to prevent accidental injuries at sharp tip.
- Do not apply extreme force to the coaxial cannula. This may cause the cannula to bend. A strongly bent coaxial cannula can impair correct and safe functionality of the product or subsequent application of compatible devices and must be replaced.
- The coaxial cannula designed for single use and must not be reused. Reuse and re-sterilization impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- Used coaxial cannula must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations.

MR safety information

MR Conditional

Non-clinical testing and simulation has demonstrated the BIP Coaxial Cannula CC-MR is MR Conditional.

A patient can be safely scanned with this product in an MR system meeting the following conditions. Failure to follow these conditions may result in injury to the patient:

Parameter	Condition of use / information
Static Magnetic Field Strength (B0)	1,5 Tesla and 3 Tesla
Static Magnetic Field Strength (B0) Orientation	Horizontal, cylindrical bore
Maximum Spatial Field Gradient (SFG)	138 T/m (13.800 G/cm)
RF Polarization	Circular polarized (CP)
RF Transmit Coil	Integrated whole body transmit RF coil
RF Receive Coil	Any receive coil may be used.
MR System (RF) Operating Modes or Constraints	Normal operating mode
B ₁ ⁺ _{rms}	Normal operating mode Note: B ₁ ⁺ _{rms} shall be used on all MR systems with this limitation parameter. Use SAR only on MR systems not providing B ₁ ⁺ _{rms} limitation.
Whole-body averaged (WBA) specific absorption rate (SAR)	2 W/kg (Normal operating mode)
Product configuration	The BIP Coaxial Cannula CC-MR can undergo an MRI scan in all configurations (with and without inner stylet).
Scan duration and wait time	15 minutes of continuous RF exposure (a sequence or back-to-back series/scan without breaks) followed by a wait time of 4 minutes before resuming scanning
MR image artifact	The presence of this product may produce an image artifact. Some manipulation of scan parameters may be needed to compensate for the artifact.

Precautions

- The warranty applies exclusively to the proper intended use of the product.
- The manufacturer does not recommend any particular coaxial cannula size. Healthcare professional users are responsible for selection of proper size depending on patient to be diagnosed.
- Please also note the operating instructions of compatible biopsy devices and compatible biopsy needles or marking clips. The manufacturer assumes no liability and warranty for use with non-compatible medical devices.

- Before use inspect each product for damages at tip or bent cannula that would prevent proper functionality. Do not use damaged or bent products.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.
Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

Environmental storage conditions

Temperature		10 – 25 °C
Humidity		20 – 50 % r.h (non-condensing)
Shelf-life		3 years

NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

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©2026 BIP GmbH, Germany, all rights reserved.