BIP Mam-O-Guide®

Localization Wire

de - Gebrauchsanweisung.....	4
en - Instructions for Use	6
fr - Mode d'emploi	9
it - Istruzioni per l'uso	11
es - Instrucciones de uso	14
pt - Instruções de utilização	16
nl - Gebruiksaanwijzing	19
da - Brugsanvisning	21
sv - Bruksanvisning	24
fi - Käyttöohje	26
cs - Návod k použití	29
sk - Návod na použitie	31
sr - Uputstvo za upotrebu	33
hu - Használati útmutató	36
ro - Instrucțiuni de utilizare	38
bg - Инструкции за употреба	41
el - Οδηγίες χρήσης	43

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ENGLISH

Intended purpose

The localization wire is intended for locating soft tissue lesions. It is to be placed preoperatively under imaging procedures, to guide the surgeon to the lesion during surgical intervention. It is intended for single-use and for short-term invasive application (maximum 24 hours).

• Indications

Indications are anomalies in breast tissue.

• Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder, anticoagulant treatment and allergy to nickel. In case of increased risks, application depends on decision of healthcare professionals.

• Residual risks and side effects

Residual risks and undesirable side effects are allergic reaction to nickel, wire dislocation and general puncture-related complications such as bleeding, infection, injury to surrounding tissue and pain.

• Intended Users

Intended users are health care professionals only (physicians).

• Intended patients

Intended patients are adults.

Clinical benefits

Universal application for all diagnostic procedures - ultrasound, X-ray (mammography) or magnetic resonance imaging MRI. Optimized anchoring in tissue with minimal risk of injury due to ring-shaped anchoring element.

For easy repositioning the localization wire can be retracted into the cannula and repositioned.

Double twisting of localization wire enables optimal reflexion and visibility to ultrasound and protects against accidental cutting of the wire.

By turning the application cannula for fine adjustment, the deployment direction of the anchoring element can be determined before the localization wire is pushed out.

BIP Mam-O-Guide®

Localization Wire

Palpable markings in front of the anchoring element enable the operator to promptly expand the incision. Optimized cannula tip facilitates penetration into tissue.

Description

The BIP Mam-O-Guide® consists of two components – an application cannula and a localization wire with anchoring element that lies in stretched form in the distal end of the application cannula.

The localization wire is made of twisted Nitinol (nickel-titanium shape memory wire) and has two short, visually and tactilely distinct indicating sections (= wider twisted area), see Fig. 4. The proximal section signals full deployment of the anchoring element and the distal section signals approach to the anchoring element during surgical intervention.

The distal end of the BIP Mam-O-Guide® forms a ring-shaped anchoring element while applied.

The application cannula has a sharp tip, centimeter depth marks for direct visual assessment of insertion depth and a small handle. The application cannula is made of medical grade stainless steel or titanium (MR-compatible variants).

The BIP Mam-O-Guide® is available in different lengths (see Tab. 1) and provided with a protective hose.

The depth stopper (see Tab. 1) can be used for pre-adjustment of penetration depth.

The BIP Mam-O-Guide® is sterile, single-packaged in peel pouches. Respective sales units are additionally packaged in protective cardboard box.

Compatibility

The BIP Mam-O-Guide® can be used with conventional imaging using ultrasound, X-ray (mammography) or magnetic resonance imaging MRI.

If the BIP Mam-O-Guide® shall be used in magnetic resonance devices (MR-systems), ensure that a MR-compatible variant is used (see Tab. 1).

Application

Preparation of the patient

The procedure is to be performed applying aseptic techniques. The lesion localization normally takes place through healthy, uninjured skin. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be used according to the discretion of the user.

Application of the medical device

Check sterile packaging and the product for expiration date and visual damages before use.

Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly 

Remove protective hose from the application cannula and place aside.

Insertion of the localization wire:

Select the site of entry for the application cannula and insert the cannula into the tissue under the control of a suitable imaging procedure up to the lesion. The depth marks facilitate exact positioning. Rotation of the application cannula around its axis during insertion can prevent unintentional bending.

The ring-shaped deployment of the anchoring element is always aligned in the opposite direction to the dot marking on the cannula handle, see Fig. 3. The alignment can be freely selected by rotating the cannula beforehand.

To anchor the localization wire, hold the application cannula firmly at the plastic handle and do not pull the cannula back. Grasp the localization wire at the proximal end of the cannula and advance slowly and gradually. The ring-shaped anchoring element is fully deployed when the proximal indicating section has just disappeared into the cannula, see Fig. 3. Fixation of anchoring element can be controlled by slight and careful jerking at the wire. After placing the localization wire, hold the wire in position and carefully remove the application cannula.

Fix the proximally protruding end of the localization wire to the breast and protect it with adequate wound care at the discretion of the user until surgical intervention.

Repositioning or removal of the localization wire:

If the localization wire is to be repositioned or removed, hold the localization wire firmly and advance the application cannula slowly and carefully until stop of the anchoring element. The localization wire can then be retracted into the cannula. The localization wire is fully retracted into the cannula when the indicating section of the localization wire is completely visible at the proximal end of the cannula, see Fig. 2. The cannula including the localization wire can then be repositioned or removed.

Disposal:

After use, dispose of the application cannula and localization wire in accordance with local provisions and safety regulations.

REACH declaration at <https://www.bipmedical.com/en/about-us/regulations/> must be observed.

Warnings

- When used on patients with breast implants, the BIP Mam-O-Guide® must not be placed in the implant and the implant must not be damaged during application.
- The BIP Mam-O-Guide® must not be used and must be disposed of if the expiration date has elapsed or the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to insertion to prevent accidental injuries at sharp tip.
- Do not apply extreme force to the application cannula. This may cause the cannula to bend. A strongly bent cannula can impair correct and safe functionality of the product and must be replaced.
- The BIP Mam-O-Guide® is intended for single use and must not be reused. Reuse and re-sterilization impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- If the tissue is removed using RF surgery, make sure that the BIP Mam-O-Guide® localization wire does not come into direct contact with the RF knife. In the event of contact, high temperatures can arise at punctiform locations and the wire could be severed as a result.
- Used product must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations.

BIP Mam-O-Guide® Localization Wire

MR safety information

MR Conditional

Non-clinical testing and simulation has demonstrated that the BIP Mam-O-Guide® MR is MR Conditional.

A patient may be safely scanned with the product in a MR system under the following conditions. Failure to follow these conditions may result in injury to the patient:

Parameter	Condition of use / information
Static Magnetic Field Strength (B0)	1,5 Tesla and 3 Tesla
Static Magnetic Field Strength (B0) Orientation	Horizontal, cylindrical bore
Maximum Spatial Field Gradient (SFG)	138 T/m (13.800 G/cm)
RF Polarization	Circular polarized (CP)
RF Transmit Coil	Integrated whole body transmit RF coil
RF Receive Coil	Any receive coil may be used.
MR System (RF) Operating Modes or Constraints	Normal operating mode
B ₁ ⁺ _{rms}	Normal operating mode Note: B ₁ ⁺ _{rms} shall be used on all MR systems with this limitation parameter. Use SAR only on MR systems not providing B ₁ ⁺ _{rms} limitation.
Whole-body averaged (WBA) specific absorption rate (SAR)	2 W/kg (Normal operating mode)
Product configuration	A localization wire can undergo an MRI scan in all configurations (with and without application cannula). The localization wire shall not be connected to any external devices.
Scan duration and wait time	15 minutes of continuous RF exposure (a sequence or back-to-back series/scan without breaks) followed by a wait time of 4 minutes before resuming scanning
MR image artifact	The presence of this product may produce an image artifact. Some manipulation of scan parameters may be needed to compensate for the artifact.

Precautions

- The warranty applies exclusively to the proper intended use of the product.
- Before use inspect each product for damages at application cannula or bent cannula that would prevent proper functionality. Do not use damaged or bent products.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.
Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

Environmental storage conditions

Temperature		10 – 25°C
Humidity		20 – 50% r.h (non-condensing)
Shelf-life		3 years

NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

Technical Data / Available variants

Shape / Dimensions	Localization wire with ring-shaped anchoring element Outer diameter localization wire: Ø 0.6 mm Diameter anchoring element: Ø 6.6 mm Length localization wire: see Tab. 1
Variants	see Tab. 1
Sterilization process	Irradiation / Gamma

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©2026 BIP GmbH, Germany, all rights reserved.