de -	Gebrauchsanweisung	3
en -	Instructions for use	20
fr -	Mode d'emploi	36
it -	Istruzioni per l'uso.....	52
es -	Instrucciones de uso.....	69
nl -	Gebruiksaanwijzing	86
pl -	Instrukcja użytkowania.....	103
cs -	Návod k použití.....	120
sk -	Návod na použitie.....	135
hu -	Használati utasítás	151
ro -	Instrucțiuni de utilizare	168
bg -	инструкция за употреба	184
el -	Οδηγίες χρήσης.....	201

BIP Biomed. Instrumente & Produkte GmbH

Am Brand 1

82299 Türkenfeld

Germany

Tel.:+49(0)8193-93180

Fax:+49(0)8193-6548

info@bipmedical.com

www.bipmedical.com

CE 0123

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ENGLISH

These instructions for use are valid for the BIP EvoCore® biopsy device with article number EC2215.

Intended purpose

The biopsy device is intended for removal of soft tissue samples for histological examination.

It is intended to be used only with a compatible disposable biopsy needle.

The biopsy device is reusable and not intended for patient contact.

• Indications

Indications are anomalies in soft tissues as breast, prostate, liver, kidney, lung or thyroid.

• Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder and anticoagulant treatment. In case of increased risks, application depends on decision of healthcare professionals.

• Residual risks and side effects

Residual risks and undesirable side effects are tumor cell seeding due to displacement of tissue or cells and general puncture-related complications such as bleeding, infection, injury to surrounding tissue/organ and pain, as well as pneumothorax and air embolism if used on lung.

• Intended users

Intended users are health care professionals only.

• Intended patients

Intended patients are adults.

Clinical benefits

Fully-automatic use for easy application.

Semi-automatic use for precise positioning of biopsy needle notch.

Adjustable penetration depth – 22 mm (standard) / 15 mm (reduced).

Enables single-handed use.

Powerful spring-loaded system allows high speed needle movement and excellent biopsy sample results.

Tissue removal without forward throw.

Patented needle protection ensures safe handling and prevents possible contamination between device and cannula.

Description

The BIP EvoCore® EC2215 is an automatic spring-operated reusable biopsy device for use in combination with compatible biopsy needles. It offers a semi-automatic and fully-automatic mode with two adjustable penetration depths (22 mm / 15 mm) each. The biopsy device has a coated aluminum housing and is designed for single-handed operation. It is non-sterile, single-packaged and must be reprocessed before use. If used as intended, the BIP EvoCore® biopsy device does not have patient contact.

Compatibility

The BIP EvoCore® biopsy device is intended for use only in combination with BIP EvoCore® biopsy needles: Identification EC[xx][yyy][z] (Example: EC14100T)

xx: Diameter in G

yyy: Length in mm

z: Version, e.g. T for WingTip-Trocar® tip

The BIP EvoCore® biopsy device and compatible needles can be used with conventional imaging using ultrasound or X-ray (mammography or computer tomography CT).

Application

The environmental temperature for application is +10 °C to +40 °C.

1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. The biopsy normally takes place through healthy, uninjured skin. When using biopsy needles with a large diameter, an incision should be made beforehand with a scalpel according to the discretion of the user. In addition to the requisite

explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

2. Preparation of the biopsy device

The biopsy device must be reprocessed and inspected for visible damage prior to use (cracks, flaking of the coating, pitting). Do not use the device if damages or defects are recognisable.

3. Cocking the biopsy device

Hold the device in your hand for this purpose. Using your index finger and middle finger, press the cocking handle (1) back into the housing twice until an audible clicking sound can be heard.

With the first stroke the cannula slide is cocked and with the second stroke the needle slide is cocked.

The current cocking state can be viewed at the indicator holes (5). If the right-hand indicator is red, the cannula slide is now cocked – the sample notch is open. If the left-hand indicator is red, the cannula slide and needle slide are cocked and the biopsy device is ready for use.

If both indicators are silver, the biopsy device is fully released. Refer also to the diagrams on the device cover (2).

4. Functional test

It is recommended to perform a functional test of the device before insertion of the biopsy needle.

5. Insertion of the biopsy needle

The device must first be cocked again, because the biopsy needle can only be inserted in the fully cocked state. Then the housing cover (2) is opened by gently pressing the button (3). The biopsy needle (4) is inserted on the two protruding guide pins, wherein it must be ensured that the biopsy needle itself is positioned on the hinge side. Then the housing cover of the biopsy device is closed again.

6. Setting the penetration depth

The standard setting is 22 mm. To select 15 mm, the slide control (6) must be pushed upwards. The penetration depth can only be switched while the device is in the fully cocked state. The slide control must be pushed completely up to the stop.

7. Removing the protective hose

Gently rotate grip of the protective hose 90° in either direction until it noticeably breaks off. Then, remove protective hose from needle and place aside. Only remove the protective hose immediately prior to performing biopsy to prevent accidental injuries at sharp tip.

8. Insertion into tissue

The insertion of the needle into the body should be carried out under the control of a suitable compatible imaging procedure.

When using the biopsy needle without a coaxial cannula, the needle is directly inserted and positioned for tissue removal.

When using the biopsy needle in combination with a coaxial cannula, the coaxial cannula must be inserted first to target and the stylet must be pulled out. The biopsy needle is then inserted through the cannula until stop and positioned for tissue removal.

The use of BIP coaxial cannula is optional, but recommended in the case of multiple biopsies in order to minimize trauma and prevent risk of tumor seeding.

For further details on biopsy needles and compatible coaxial cannulas refer to instructions for use of BIP EvoCore® biopsy needles.

9. Performing the biopsy

After the successful positioning of the BIP EvoCore® biopsy needle and without interim change in position, the shot is triggered, for which purpose the user can choose between semi-automatic and fully-automatic mode:

Semi-automatic mode

Press the **Single** button (7), wherein only the inner needle is shot. The outer cannula is shot with the actuation of the **Automatic** button (8). This mode can significantly increase the biopsy yield, depending on the tissue characteristics.

The biopsy needle must never be pulled out of the body without pressing the **Automatic** button (8), because a tissue sample could not be removed and there would be the risk of tissue transfer into the punch channel.

Fully-automatic mode

Press the **Automatic** button (8), wherein the inner needle and outer cannula are triggered in very short succession. The biopsy process is completely performed with a push of the **Automatic** button.

10. Removal of the tissue sample

For the removal of the tissue sample, the outer cannula is retracted by cocking once with the cocking handle (1) and the tissue sample can then be extracted.

Therefore, it is not necessary to remove the biopsy needle from the device for the removal of the tissue sample.

If additional biopsies should take place after the removal of the tissue sample, the device is ready for use after re-cocking

11. Removal of the biopsy needle

If no further tissue samples are to be taken after the removal of the tissue sample, the cocking handle (1) must then be pulled back. Then the housing cover can be opened by gently pressing the button (3) and the biopsy needle can be removed. The biopsy needle must be disposed of in accordance with local provisions and safety regulations.

After removal of the biopsy needle, press the **Automatic** button (8) to switch the device to a fully released state (both indicator holes are silver).

Warnings

- The biopsy device must be prepared in accordance with these instructions before each use
- The length of the biopsy needle must be selected so that the biopsy device does not come into contact with the patient.
- Biopsy device and needles are not MR compatible. For use in MRI environment ensure that the biopsy procedure is only applied outside 5 Gauss line and MR compatible coaxial cannula is used for determination of target position.
- Do not make "empty shots" with inserted biopsy needle. Activating the biopsy device in the air, not in tissue, can result in damage to the cannula tip and can impair the correct and safe performance of the biopsy needle.
- The cocking handle must not be pulled during the release (shot) to ensure the correct and safe performance.

Proper cleaning and disinfection

General principles

Before each use of the biopsy device, BIP GmbH recommends a complete reprocessing of the biopsy device. Simplified cleaning/disinfecting of the biopsy device are the user's responsibility and take place according to the user's discretion. This applies in particular when using the device for the first time after delivery, because all biopsy devices are delivered unsterile.

The concentrations, temperatures and exposure times specified by the manufacturer of the cleaning and disinfecting agents as well as the requirements for rinsing must be strictly observed. Use only freshly prepared solutions, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. highly purified water HPV) and only medical compressed air (acc. to European Pharmacopoeia) for drying.

Please observe, in the scope of your responsibility, for the sterility of the biopsy device to be used,

- that only sufficiently validated product-specific methods for cleaning/disinfecting and sterilisation are used as a fundamental rule.
- that the used equipment (disinfector, steriliser) is regularly maintained and inspected.
- that the validated parameters are observed for each cycle.

Please also observe the valid legal regulations in your country as well as the hygienic regulations of the medical practice and/or hospital. This particularly applies for the various specifications with respect to an effective inactivation of prions.

If possible, avoid extensive soiling of the biopsy device during and after use.

Use of the biopsy device after prior cleaning and/or disinfecting without subsequent sterilisation takes place under the sole responsibility of the user.

Cleaning and disinfection

Basics

If possible, a machine process (disinfector) should be used for the cleaning and disinfecting. A manual process should only be used if a machine process is not available due to the significantly lower effectiveness and reproducibility.

Pre-treatment must be performed in both cases.

Conditioning instructions as per EN ISO 17664-1

Warning notices:	do not use severely basic or acid-containing cleaners no concentrated chlorhexidine gluconate or other media containing alcohol.
Limit on re-conditioning:	none

Risk assessment and evaluation of medical products prior to conditioning

The type and scope of conditioning are related to the application of the medical product. For this reason, the operator shall be responsible for the correct assessment of medical products and for specifying the type and scope of the conditioning (see KRINKO/BfArM recommendation*, item 1.2.1 risk evaluation and assessment of medical products prior to conditioning). Based on this user-dependent assessment, the operator may specify which of the conditioning processes listed in this conditioning instruction must be completed.

* Bundesgesundheitsblatt 2012 55:1244-1310; DOI 10.1007/s00103-012-1548-6.

Pretreatment

For the selection of the cleanser and disinfectant to be used, it must be ensured

- that a disinfectant with certified validity (e.g. VAH/DGHM or FDA approval and/or CE marking) is used and that it is compatible with the cleanser which is used and
- that the chemicals which are used are compatible with the biopsy device (see chapter "Material Resistance").

Please note that the disinfectant used in the pretreatment is only used for personal protection and cannot replace the disinfection step to be carried out later after cleaning.

Instructions

Location of use	Remove and dispose of the biopsy needle.
Storage and transport	Transport the contaminated device in a closed container. We recommend preparing the device again within 2 hours after use.
Conditioning	Personal protective equipment (wear gloves, water-repellent protective coat, face protection mask or safety glasses, follow the instructions of the disinfectant and cleaner manufacturer.
Pretreatment	<ol style="list-style-type: none">1) Set the device to the 22 mm position with the cover open, since this enables the best possible access to the interior for cleaning.2) Submerge the device completely in tap water for at least 5 minutes.3) Flush out the device for at least 5 minutes under running tap water while moving the mechanism while tensioning and releasing the device 5 times.
Manual cleaning	Equipment: Cleaner (multi-stage enzymatic cleaner featuring manufacturer's permission for anodised aluminium), soft brush, tray for cleaner. <ol style="list-style-type: none">1) Apply the cleaner according to the manufacturer's specifications.2) Submerge the pre-set device with the flap open and the cover down completely in the cleaning solution.

	<p>3) Allow to work for at least 15 minutes.</p> <p>4) During the reaction time, close the flap of the device, tension the device, release it, and open the cover. Repeat all steps at least 4 times.</p> <p>5) Afterwards, brush the device for at least 1 minute while tensioning the white tensioning slider so that the brush reaches all accessible areas.</p> <p>6) Remove the device from the cleaning solution.</p> <p>7) Rinse the device with the flap opened for at least 1 minute under running tap water.</p> <p>8) Check the cleanliness of the device and repeat the previously indicated steps if fouling is visible.</p>
Manual disinfection	<p>Equipment: Peroxide-based disinfectant (peracetic acid), disinfectant tray, medical compressed air (acc. to European Pharmacopoeia).</p> <p>1) Apply the disinfectant solution according to the manufacturer's specifications.</p> <p>2) Submerge the opened device with the cover down completely in the disinfectant solution.</p> <p>3) During the reaction time close the cover while submerged, tension the device, release it, and open the cover. Repeat all steps at least 4 times.</p> <p>4) Allow the disinfectant to work according to the manufacturer's specifications.</p> <p>5) Remove the device from the disinfectant solution while open.</p> <p>6) Flush out the device for at least 1 minute under running tap water while tensioning and releasing the device 5 times.</p> <p>7) Wipe the device with a lint-free disposable cloth and dry with medical compressed air.</p>
Mechanical cleaning and disinfecting	<p>Complete pretreatment prior to mechanical cleaning and thermal disinfection (see pretreatment).</p> <p>Equipment: Washer-disinfector (WD) as per EN ISO 15883-1/ -2 with thermal program (temperature 90 °C to 95 °C); mild-alkali cleaner</p>

with manufacturer's permission for anodised aluminium, medical compressed air (acc. to European Pharmacopoeia).

- 1) Set the device to the 22 mm position with the cover open and insert it in the washer-disinfector with the opening facing downwards.
- 2) Close the washer-disinfector and start the thermal program, see the following table for the program process.
- 3) After the program completes, remove the device.
- 4) Check if the device is completely dry, wipe it off with a lint-free disposable cloth as required, and dry it with compressed air.

Prog. step	Time	Water	Dosage	Temperature
Pre-rinse	5 min	cold		--
Dose Cleaner			As per manufacturer's specifications	As per manufacturer's specifications
Clean	10 min	demineralised		55 °C
Rinse	2 min	demineralised		--
Disinfection	5 min	demineralised	As per manufacturer's specifications	90 °C
	95 sec			95 °C
Dry	15 min			90 °C - max. 120 °C

Maintenance, control, and inspection

Complete visual inspection for damage, corrosion, and wear on the device. Check the functionality of the triggering mechanism. Lubricate all moving parts of the device with a silicone-free, steam-sterilised lubricant with medical approval according to the instructions of the manufacturer.

Packing	<p>Packing in sterile barrier system as per EN ISO 11607 or DIN 58953</p> <p>1) The device must be in position 22 mm and <u>must not be cocked</u>. Both indicators on the housing cover must be silver.</p> <p>2) Pack the device.</p>
Sterilisation	<p>-Steriliser as per EN 13060 and/or EN 285.</p> <p>-Sterilisation with moist heat according to EN ISO 17665-1/-2.</p> <p>-Fractionated vacuum process, 134 °C, 200 – 2000 mbar, hold time of at least 3 minutes.</p> <p>1) Place the packed device in the sterilisation chamber.</p> <p>2) Start the program.</p> <p>3) After the program completes, remove the device and allow it to cool off.</p> <p>4) Check that the packing is dry and free of damage.</p>
Storage	Storage and storage duration according to the specifications of the conditioner.
Additional information	none

The instructions indicated above were validated by the manufacturer of the medical product as a suitable process for re-conditioning the medical product for re-use. It remains the responsibility of the user to ensure that the currently implemented conditioning process of the device achieves the desired and necessary results corresponding with use of equipment, materials, and personnel. This requires the process to be validated and checked routinely. Accordingly, any deviation on behalf of the user from the indicated instructions with regard to effectiveness and potentially unfavourable consequences must be evaluated and assessed.

Storage

After sterilisation, the biopsy device must be stored in the sterilisation packaging at room temperature in a dry, dust-free location protected from light.

Material resistance

When selecting the cleanser and disinfectant, make sure that they do not contain the following components:

- organic, mineral and oxidising acids (minimum permissible pH value 5.5)
- stronger bases (maximum permissible pH value 9.5, neutral/enzymatic cleaner recommended)
- organic solvents (e.g. alcohols, ether, ketone, benzene)
- oxidants (e.g. hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

In case an ultrasound bath is used, the anodised coating could bleach and the housing labels could erode. However, this does not affect performance.

Never clean the biopsy device with metal brushes or steel wool.

Re-usability

The biopsy device is not subject to any maximum validated conditioning cycles. The requirement for safe functioning is correct completion of conditioning and regular care of the device. The manufacturer recommends having the biopsy device maintained by the manufacturer after 100 conditioning cycles.

Instructions for lubrication and care

The factory lubricant can be removed from by the cleanser over the course of time. In order to ensure the performance and service life of the device, the biopsy device must be oiled at the slide guides (A) and (B) and the springs opposite, both in the cocked and uncocked state each time it is prepared. Lubricating is done with a silicone-free lubricant compatible with steam sterilization. The manufacturer's instructions for the respective lubricant must be followed. Aesculap Sterilit oil spray JG600 has been proven to be effective.

Troubleshooting

Problem	Solution	Alternative
The device cannot be cocked or only with difficulty.	<ul style="list-style-type: none"> ◇ Hold the device in your hand according to the specifications for the general description of use. ◇ Pull the cocking handle back until a clear click can be heard. ◇ Make sure that you do not press a trigger button at the same time. ◇ Make sure that no foreign particles are blocking the cocking mechanism in the interior of the device. 	Send the device to BIP GmbH for inspection.
The cover does not open.	<ul style="list-style-type: none"> ◇ Make sure that the cover is not canted and/or the opening is not blocked. ◇ Press the button to the appropriate depth for the cover unlocking (3). ◇ Clean and lubricate the device. 	Send the device to BIP GmbH for inspection.
The biopsy needle can not be inserted.	<ul style="list-style-type: none"> ◇ Make sure that the device is fully cocked. ◇ Make sure that the holders for the positioning of the biopsy needle are not deformed. ◇ Make sure that only a BIP EvoCore® biopsy needle is used. 	Send the device to BIP GmbH for inspection.
The cover does not close.	<ul style="list-style-type: none"> ◇ Check whether the cover can be closed without the biopsy needle inserted. ◇ Make sure that the cover is not canted and/or the closing is not blocked. 	Send the device to BIP GmbH for inspection.

	<ul style="list-style-type: none"> ◇ Make sure that the BIP EvoCore® biopsy needle is positioned correctly. ◇ Clean and lubricate the device. 	
Biopsies cannot be performed in Automatic mode and/or only with difficulty.	<ul style="list-style-type: none"> ◇ Make sure that the device is fully cocked. ◇ Make sure that the slide control is engaged at the respective stop (22 mm or 15 mm). ◇ Make sure that you completely press the automatic button. ◇ Make sure that no foreign particles are blocking the mechanism in the interior of the device. ◇ Clean and lubricate the device. 	Send the device to BIP GmbH for inspection.
Biopsies cannot be performed in Single mode and/or only with difficulty.	<ul style="list-style-type: none"> ◇ Make sure that the device is fully cocked. ◇ Make sure that the slide control is engaged at the respective stop (22 mm or 15 mm). ◇ Make sure that you completely press the trigger buttons. ◇ Make sure that no foreign particles are blocking the mechanism in the interior of the device. ◇ Clean and lubricate the device. 	Send the device to BIP GmbH for inspection.

<p>There is no tissue sample in the sample notch.</p>	<ul style="list-style-type: none"> ◇ Make sure that you press the Automatic button after pressing the Single button. ◇ Check (without biopsy needle inserted) whether both cannula slides are activated in the corresponding modes. In this connection, pay attention to the status display on the indicator holes. ◇ Make sure that the slide control is engaged at the respective stop (22 mm or 15 mm). 	<p>Send the device to BIP GmbH for inspection.</p>
<p>The penetration depth cannot be adjusted.</p>	<ul style="list-style-type: none"> ◇ Make sure that the device is fully cocked. ◇ Clean and lubricate the device. 	<p>Send the device to BIP GmbH for inspection.</p>

Maintenance and repair

Maintain the cleaned and disinfected device after 100 cleaning cycles.

	<p>Maintenance and repair must only be performed by the manufacturer or by qualified personnel expressly authorised by the manufacturer.</p>
--	---

Observe the following periods:

Period	Relevant parts	Person performing tasks
<p>Maintenance → every 100 cleaning cycles</p>	<p>Defined safety-relevant wear parts</p>	<p>Manufacturer or qualified personnel expressly authorised by the manufacturer</p>

Warranty and precautions

- BIP guarantees that the BIP EvoCore® biopsy device will remain free of defects when used as intended for a period of two years.
The submission of a receipt of purchase indicating the seller and date of purchase is required for warranty claims.
- BIP assumes no liability and warranty for the following:
 - Non-observance of the instructions for use
 - Operating errors
 - Improper use or improper handling
 - Use with non-compatible medical devices
 - Third-party interference by unauthorised persons into the instrument for repair purposes
 - Transport damage due to improper packaging for returns
 - Neglected maintenance
 - Ageing due to operation and normal wear.This includes the following components, for example:
 - Springs
 - Articles intended for single use, etc.
 - Use of non-original spare parts
- BIP is not liable for consequential damages from defects, insofar as they due to intent or gross negligence or for injury to life or limb due to ordinary negligence.
- BIP reserves the right, according to its own discretion, to rectify defects, to deliver a fault-free item or to appropriately reduce the purchase price.
- With the rejection of a warranty claim, the costs for transport and return transport are not assumed.
- Please also note the operating instructions of compatible biopsy needles and accessories.
- According to the European Medical Device Regulation (2017/745) BIP is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to BIP and the competent authority of the member state in which the user and/or patient is established.
Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

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©2023 BIP GmbH, Germany, all rights reserved.