# BIP O-Twist-Marker Marking Clip


de	-	Gebrauchsanweisung	4
en	-	Instructions for use	7
fr	-	Mode d'emploi	9
it	-	Istruzioni per l'uso	12
es	-	Instrucciones de uso	14
nl	-	Gebruiksaanwijzing	17
pt	-	Instruções de utilização	19
da	-	Brugsanvisning	22
sv	-	Bruksanvisning	25
fi	-	Käyttöohje	27
cs	-	Návod k použití	30
sk	-	Návod na použitie	32
hu	-	Használati útmutató	35
ro	-	Instrucțiuni de utilizare	37
bg	-	Инструкции за употреба	40
el	-	Οδηγίες χρήσης	43

BIP Biomed. Instrumente & Produkte GmbH  
Am Brand 1  
82299 Türkenfeld  
Germany  
Tel.: +49(0)8193-93180  
Fax: +49(0)8193-6548  
[info@bipmedical.com](mailto:info@bipmedical.com)  
[www.bipmedical.com](http://www.bipmedical.com)

CE 0123

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## ENGLISH

### Intended purpose

The implantable clip is intended for marking of soft tissue for monitoring of lesions after biopsies, before or during chemotherapy and in radiotherapy treatment planning. It is intended for single-use and for long-term application.

#### • Indications

Indications are anomalies in breast tissue including axillary lymph nodes.

#### • Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder, anticoagulant treatment and allergy to nickel.

In case of increased risks, application depends on decision of healthcare professionals.

#### • Residual risks and side effects

Residual risks and undesirable side effects are allergic reaction to nickel, marker dislocation and general puncture-related complications such as bleeding, infection, lymphedema, injury to surrounding tissue and pain.

#### • Intended users

Intended users are health care professionals only.

#### • Intended patients

Intended patients are adults.

### Clinical benefits

Universal application for all diagnostic procedures - ultrasound, X-ray (mammography, computer tomography CT) or magnetic resonance imaging MRI.

Twisted, ring-shaped clip element enables favorable reflection and visibility in ultrasound.

Secure anchorage in the tissue with minimal dislocation.

Universal use following conventional core or vacuum biopsy procedures, independent of biopsy system.

Applicator with sharp cannula for direct percutaneous application of marker.

### Description

The BIP O-Twist-Marker consists of two components – an applicator and a marking clip that lies in stretched form in the distal end of the application cannula. The clip is made of braided Nitinol (nickel-titanium shape memory alloy) and forms a ring when placed in tissue. The applicator consists of an outer cannula with sharp tip and centimeter depth marks for direct visual assessment of insertion depth (except type SRX), and an inner movable mandrel for advancing the marking clip, both connected to plastic handles. The cannula is made of medical grade stainless steel or titanium (MR-compatible variant), the inner mandrel is made of nickel-titanium alloy. The handles of the application cannula and the mandrel are secured against unintentional actuation by a release protector.

The applicator is available in different lengths (see Tab. 1) and is provided with a protective hose.

Type S / SA is the straight version, which is intended for direct advancement through the tissue.

Type R / RA / SRX is provided with a radius at the distal end and is intended for application through the existing sheath of a vacuum biopsy system.

The depth stopper (see Tab. 1) can be used for pre-adjustment of application depth.

Type SRX contains additionally separately packed adapter SRX for compatible biopsy device (see Fig. 3).

The BIP O-Twist-Marker is sterile, single-packaged in peel pouches. Respective sales units are additionally packaged in protective cardboard box.

### Compatibility

The BIP O-Twist-Marker types S, S-MR, SA can be used in combination with compatible BIP Coaxial Cannulas (see Tab. 1).

The BIP O-Twist-Marker type SRX is intended for use only in combination with BARD® EnCor Enspire™ Biopsy System 7 G / 10 G.

The BIP O-Twist-Marker can be used with conventional imaging using ultrasound or X-ray (mammography, computer tomography CT) or magnetic resonance imaging MRI.

If the BIP O-Twist-Marker shall be used in magnetic resonance devices (MR-systems), ensure that the MR compatible variant OTM3.0S137MR is used (see Tab. 1).

### Application

#### Preparation of the patient

The procedure is to be performed applying aseptic techniques. The tissue marking normally takes place through healthy, uninjured skin. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

#### Application of the medical device

##### Type S / S-MR / R / SA / RA:

1. Open peel pouch and aseptically remove applicator along with the included marking clip.
2. Remove protective hose from application cannula and place aside.
3. **OTM3.0S / S-MR / OTM3.0SA:** Under the control of suitable imaging procedure, advance the applicator directly or through compatible coaxial cannula into the tissue up to the intended marking point, rotating it at the same time.  
NOTE: With OTM3.0S the depth stop must be removed before use in combination with a coaxial cannula.  
**OTM3.0R / OTM3.0RA:** Under the control of suitable imaging procedure, advance the applicator through the sheath of a vacuum biopsy system up to the intended marking point. Make sure that in sheaths with a side window the tip of the cannula is at the same angle of rotation to the outlet window. NOTE: The bending direction of cannula as well as deployment orientation of marking clip is indicated by the marking or the blue arrow on the handle (see "Marking" Fig. 4). To prevent migration of the marker when positioned after a vacuum biopsy, the marker should not be deposited directly into the biopsy cavity but instead, it should be fixed by inserting the bent cannula tip into the tissue bordering the biopsy cavity.
4. Remove release protector between cannula and mandrel (see Fig. 4).
5. Slowly advance the mandrel until stop, applying uniform pressure on the proximal end of the handle.
6. Remove the applicator and dispose of in accordance with local provisions and safety regulations.

# BIP O-Twist-Marker Marking Clip

## Type SRX

1. A precondition is that the BARD® EnCor Enspire™ biopsy needle must be correctly positioned with an open sample notch.
2. Remove the closure on the tissue collection container of the BARD® EnCor Enspire™ biopsy needle.
3. Open peel pouch and aseptically remove the adapter SRX ①.
4. Insert the adapter SRX into the opening at the proximal end of the BARD® EnCor Enspire™ biopsy needle until stop. The blue seal must no longer be visible. The correct position of the adapter SRX is reached when a clearly noticeable resistance has been overcome.
5. Open peel pouch and aseptically remove applicator along with the included marking clip ②.
6. Remove protective hose from application cannula and place aside.
7. Completely insert the applicator into the proximal end of the adapter SRX until stop. In doing so make sure that the red mark on the BARD® EnCor Enspire™ biopsy needle and the blue arrow on the application handle (see "Marking" Fig. 4) line up. To prevent migration of the marker when positioned after a vacuum biopsy, the marker should not be deposited directly into the biopsy cavity but instead, it should be fixed by inserting the bent cannula tip into the tissue bordering the biopsy cavity.
8. Remove release protector between cannula and mandrel (see Fig. 4).
9. Slowly advance the mandrel until stop, applying uniform pressure on the proximal end of the handle.
10. Slightly withdraw and rotate the applicator, then rotate the BARD® EnCor Enspire™ biopsy needle 180° and close the sample notch.
11. Pull the BARD® EnCor Enspire™ biopsy needle with the applicator as a unit out of the puncture channel. Dispose of applicator in accordance with local provisions and safety regulations.

## Warnings

- The marking clip cannot be absorbed and remains permanently implanted if removal via biopsy or operation is not planned.
- When used on patients with breast implants, the BIP O-Twist-Marker must not be placed in the implant and the implant must not be damaged during application.
- The BIP O-Twist-Marker must not be used and must be disposed of if the expiration date has elapsed or the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to insertion to prevent accidental injuries at sharp tip.
- Do not apply extreme force to the applicator. This may cause the cannula to bend. A strongly bent cannula can impair correct and safe functionality of the product and must be replaced.
- Both the applicator and the marking clip are designed for single use and must not be reused. Reuse or reloading impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- Used applicator must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations.

## MR safety information

MR Conditional

### Type S-MR:

Non-clinical testing and simulation has demonstrated the application cannula of type S-MR is MR Conditional.

A patient can be safely scanned with this product in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 13.800 G/cm (138 T/m)
- Maximum force product of 225.000.000 G<sup>2</sup>/cm (225 T<sup>2</sup>/m)
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the BIP application cannula is expected to produce a maximum temperature rise of less than

- 1,6 °C (2 W/kg, 1.5 Tesla / 64 MHz) RF-related temperature increase
- 7,4 °C (2 W/kg, 3 Tesla / 128 MHz) RF-related temperature increase

after 15 minutes of continuous scanning.

Application cannulas of all other variants are not MR compatible and must only be used in MR environment outside 5 Gauss line!

### Marking clip: (all types)

A patient can be safely scanned with the marking clip in an MR system meeting the following conditions:

According to ASTM F2182 a RF-related temperature increase of < 2° C during 1 hour of continuous scanning at 1.5 Tesla / 64 MHz or 3 Tesla / 128 MHz frequencies can be assumed, since the marking clip has an overall dimension of < 2 cm. This applies to the condition that multiple marking clips are not implanted within 3 cm area.

The expected RF-related temperature increase was confirmed for a marking clip via non-clinical tests for the following scanning conditions: 0.1 ± 1.5 °C (4 W/kg, 3 Tesla / 128 MHz) after 15 minutes of continuous scanning.

## Precautions

- The warranty applies exclusively to the proper intended use of the product.
- Please also note the operating instructions of compatible products.
- Before use inspect each product for damages at applicator or bent cannula that would prevent proper functionality. Do not use damaged or bent products.
- According to the European Medical Device Regulation (2017/745) BIP is obligated to inform you of the following:  
All serious incidents that have occurred in relation to the medical device must be reported to BIP and the competent authority of the member state in which the user and/or patient is established.

Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

# BIP O-Twist-Marker Marking Clip

## Environmental storage conditions

Temperature 	10 – 25 °C
Humidity 	20 – 50 % r.h. (non-condensing)
Shelf-life 	3 years

## NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

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©2022 BIP GmbH, Germany, all rights reserved.