# BIP EvoCore<sup>®</sup> Biopsy Needle


de	-	Gebrauchsanweisung.....	4
en	-	Instructions for use .....	6
fr	-	Mode d'emploi .....	8
it	-	Istruzioni per l'uso.....	10
es	-	Instrucciones de uso.....	12
nl	-	Gebruiksaanwijzing .....	14
no	-	Bruksanvisning .....	16
pl	-	Instrukcja użytkowania .....	18
cs	-	Návod k použití.....	20
sk	-	Návod na použitie.....	22
hu	-	Használati útmutató.....	24
ro	-	Instrucțiuni de utilizare.....	26
bg	-	Инструкции за употреба .....	28
el	-	Οδηγίες χρήσης.....	30

BIP Biomedizinische Instrumente und Produkte GmbH  
Am Brand 1  
82299 Türkenfeld  
Germany  
Tel.: +49(0)8193-93180  
Fax: +49(0)8193-6548  
[info@bipmedical.com](mailto:info@bipmedical.com)  
[www.bipmedical.com](http://www.bipmedical.com)

**CE** 0123

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## ENGLISH

### Intended purpose

The biopsy needle is intended for removal of soft tissue samples for histological examination. It is intended to be used only with a compatible reusable biopsy device. It is intended for single-use and for transient invasive application.

#### • Indications

Indications are anomalies in soft tissues as breast, prostate, liver, kidney, lung or thyroid.

#### • Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system. It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder and anticoagulant treatment. In case of increased risks, application depends on decision of healthcare professionals.

#### • Residual risks and side effects

Residual risks and undesirable side effects are tumor cell seeding due to displacement of tissue or cells and general puncture-related complications such as bleeding, infection, injury to surrounding tissue/organ and pain, as well as pneumothorax and air embolism if used on lung.

#### • Intended users

Intended users are health care professionals only (physicians).

#### • Intended patients

Intended patients are adults.

### Clinical benefits

Visibility to ultrasound or X-ray (mammography, computer tomography CT).

3D rotation loop of the cannulas (special grinding at tip) for optimal biopsy quality.

WingTip-Trocar<sup>®</sup> tip (type EC-T) for optimized tissue penetration and increased target precision.

Patented needle protection ensures safe handling and prevents possible contamination between device and cannula.

### Description

The BIP EvoCore<sup>®</sup> biopsy needle is made of medical-grade stainless steel and consists of an inner stylet with 19 mm sample notch, available with beveled tip (type EC) or with WingTip-Trocar<sup>®</sup> tip (type EC-T), an outer cannula and a protective hose connected to the plastic needle holder (see Fig. 2).

The outer cannula is equipped with centimeter depth marks for direct visual assessment of insertion depth.

The BIP EvoCore<sup>®</sup> biopsy needle is sterile, single-packaged in peel pouches and is available in various diameters and lengths (see Fig. 8, Tab. 1). Sales units are additionally packaged in protective cardboard box.

### Compatibility

The BIP EvoCore<sup>®</sup> biopsy needle is intended for use in combination with BIP EvoCore<sup>®</sup> and BARD<sup>®</sup> Magnum<sup>®</sup> biopsy devices.

The BIP EvoCore<sup>®</sup> biopsy needle can be used in combination with compatible BIP Coaxial Cannulas (see Tab. 1).

The BIP EvoCore<sup>®</sup> biopsy needle can be used with conventional imaging using ultrasound or X-ray (mammography, computer tomography CT).

### Application

#### 1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. The biopsy normally takes place through healthy, uninjured skin. When using biopsy needles with a large diameter, an incision should be made beforehand with a scalpel according to the discretion of the user. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

#### 2. Cocking of the biopsy device - Fig. 3

Cock the biopsy device completely with two strokes.

#### 3. Remove the biopsy needle from sterile package

Check sterile packaging and the product for expiration date and visual damages before use.

Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly  and remove the product aseptically.

#### 4. Insertion of the biopsy needle – Fig. 4

Open housing of biopsy device and insert the needle with needle holder onto the two protruding guide pins. Close housing cover again.

#### 5. Remove protective hose - Fig. 5 to 7

Gently rotate grip of the protective hose 90° in either direction until it noticeably breaks off. Remove protective hose from needle and place aside.

#### 6. Insertion into tissue

The insertion of the needle into the body should be carried out under the control of a suitable compatible imaging procedure.

When using the biopsy needle without a coaxial cannula, the needle is directly inserted and positioned for tissue removal.

When using the biopsy needle in combination with a coaxial cannula, the coaxial cannula must be inserted first to target and the stylet must be pulled out. The biopsy needle is then inserted through the cannula until stop and positioned for tissue removal.

The use of BIP coaxial cannula is optional but recommended in the case of multiple biopsies in order to minimize trauma and risk of tumor seeding. Make sure that correct coaxial cannula (size and diameter) is used for applied biopsy needle (see Tab. 1).

#### 7. Performing the biopsy and removal of the tissue sample

Refer to the instructions for use of the respective biopsy device being used.

## 8. Disposal

Dispose of in accordance with local provisions and safety regulations.

### Warnings

- For those organs as listed in the indications, safety and performance is demonstrated on sufficient clinical data. If the device is to be used on other similar soft tissues, except those body regions listed as contraindications, safe application depends on decision of healthcare professionals.
- The biopsy needle must not be used and must be disposed of if the product is damaged or if the expiration date has elapsed or if the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to performing biopsy to prevent accidental injuries at sharp tip.
- Do not make "empty shots". Activating the biopsy device in the air, not in tissue, can result in damage to the cannula tip and can impair the correct and safe functionality of the biopsy needle.
- Do not apply extreme force to the biopsy needle, especially while stylet is extended from cannula after "single-shot" of biopsy device. This may cause the stylet or cannula to bend. A strongly bent biopsy needle can impair correct and safe functionality of the product and must be replaced.
- Biopsy needles are not MR compatible. For use in MRI environment ensure that the biopsy is only performed outside the 5 Gauss line and MR compatible coaxial cannula (see Tab. 1, CC[xx][yyy]MR) is used for determination of target position.
- The biopsy needle is designed for single use and must not be reused. Reuse and re-sterilization impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- Used biopsy needle must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations.

### Precautions

- The warranty applies exclusively to the proper intended use of the product.
- The manufacturer does not recommend any particular biopsy needle size. Healthcare professional users are responsible for selection of proper size depending on patient to be diagnosed.
- Please also note the operating instructions of compatible biopsy devices and coaxial cannulas. The manufacturer assumes no liability and warranty for use with non-compatible medical devices.
- Before use inspect each needle for damages at tip or bent cannula that would prevent proper functionality. Do not use damaged or bent needles.
- Use of biopsy needle on dense tissue (e.g. prostate) may cause stylet to bend, especially true for smaller diameter sizes and after performing multiple biopsies. Strongly bent biopsy needles must be replaced by new ones.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:  
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.  
Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

### Environmental storage conditions

Temperature		10 – 25 °C
Humidity		20 – 50 % r.h (non-condensing)
Shelf-life		3 years

#### NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

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©2023 BIP GmbH, Germany, all rights reserved.