BIP Fine Biopsy Needle

de	-	Gebrauchsanweisung		3
en	-	Instructions for use		5
fr	-	Mode d'emploi		7
it	-	Istruzioni per l'uso		9

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ENGLISH

Intended purpose

The fine biopsy needle is intended for removal of cell samples from soft tissue for cytological examination. It is intended for single-use and for transient invasive application.

• Indications

Indications are anomalies in soft tissues as prostate.

• Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder and anticoagulant treatment. In case of increased risks, the application depends on decision of healthcare professional.

• Residual risks and side effects

Residual risks and undesirable side effects are general puncture-related complications such as bleeding, infection, injury to surrounding tissue/organ and pain.

• Intended Users

Intended users are health care professionals only (physicians).

• Intended patients

Intended patients are adults.

Clinical benefits

Diagnostic extraction of cell samples

Description

The BIP Fine Biopsy Needle is made of medical-grade stainless steel and consists of two parts, an inner stylet with a central trocar tip and plastic handle and an outer cannula with sharp tip and handle piece.

The outer cannula is equipped with centimeter depth marks for direct visual assessment of insertion depth and with a female Luer-Lock connector that enables connection of syringes to the proximal end of the cannula for aspiration biopsy (see Fig. 2, 3).

The BIP Fine Biopsy Needle is sterile single-packaged in peel pouches. Sales units are additionally packaged in protective cardboard box.

Compatibility

Sterile disposable syringes with male Luer-Lock connector.

The BIP Fine Biopsy Needle can be used with conventional imaging using ultrasound or X-ray (computer tomography CT).

Application

1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. The application normally takes place through healthy, uninjured skin. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

2. Remove the fine biopsy needle from sterile package

Check sterile packaging and the product for expiration date and visual damages before use.

Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly  and remove the product aseptically.

3. Remove the protective hose

Remove protective hose from fine biopsy needle and place aside.

4. Insertion into tissue

Insert the fine biopsy needle under the control of suitable compatible imaging procedure until the stylet tip is directly in front of the area to be treated. Hold the plastic handle of the cannula and pull out the inner stylet at the handle with a slight twisting motion. The alignment of the needle tip can be checked using the raised dot on the cannula handle.

5. Performing the biopsy

Connect a syringe with a Luer-Lock connector to the handle of the cannula.

By carefully pulling back the syringe plunger manually, cells are sucked into the cannula lumen by negative pressure. After cell collection, first release the syringe plunger to reduce the negative pressure at the tip of the cannula.

Only then the fine biopsy needle should be withdrawn from the tissue, as otherwise there is a risk of aspirating other tissue cells.

Remove the syringe from the fine biopsy needle, pull back the plunger again, and fill the syringe with air. Then reconnect the syringe with the Luer-Lock connector to the handle of the cannula and apply the aspirate to a microscope slide by pushing the syringe plunger forward. If several samples are to be taken, reinsert the stylet into the cannula and reposition the fine biopsy needle using imaging techniques.

Further cell samples are taken as described above.

6. Disposal

After use, dispose of in accordance with local provisions and safety regulations.

REACH declaration at <https://www.bipmedical.com/en/about-us/regulations/> must be observed.

Warnings

- For intended purpose and those organs as listed in the indications, safety and performance is demonstrated on sufficient clinical data. If the device is to be used alternatively on other similar soft tissues, except those body regions listed as contraindications, safe application depends on decision of healthcare professionals.
- The fine biopsy needle must not be used and must be disposed of if the product is damaged or if the expiration date has elapsed or if the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to insertion to prevent accidental injuries at sharp tip.
- Do not apply extreme force to the fine biopsy needle. This may cause the cannula to bend. A strongly bent fine biopsy needle can impair correct and safe functionality of the product and must be replaced.
- Fine biopsy needles are not MR compatible. For use in MRI environment ensure that the biopsy is only performed outside the 5 Gauss line.
- The fine biopsy needle is designed for single use and must not be reused. Reuse and re-sterilization impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- Used fine biopsy needle must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations.

Precautions

- The warranty applies exclusively to the proper intended use of the product.
- Before use inspect each product for damages at tip or bent cannula that would prevent proper functionality. Do not use damaged or bent products.
- Use of fine biopsy needle on dense tissue may cause fine biopsy needle to bend, especially true after performing multiple biopsies. Strongly bent fine biopsy needles must be replaced by new ones.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

Environmental storage conditions

Temperature 	10 – 25°C
Humidity 	20 – 50% r.h (non-condensing)
Shelf-life 	3 years

NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

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©2026 BIP GmbH, Germany, all rights reserved.