Needle Guide Bushing STBG-CY

de	-	Gebrauchsanweisung	3
en	-	Instructions for use	5

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ENGLISH

Intended purpose

The needle guide bushing is intended for guiding compatible needles and cannulas.
It is intended to be used only with a compatible stereotactic system.
The needle guide bushing is intended for single-use and has no patient contact.

• Indications

Indications are defined by the compatible stereotactic system according to its application and are not given by the needle guide bushing itself.

• Contraindications

Contraindications are defined by the compatible stereotactic system according to its application and are not caused by the needle guide bushing itself.

• Residual risks and side effects

Residual risks and undesirable side effects are defined by the compatible stereotactic system according to its application and are not caused by the needle guide bushing itself.

• Intended Users

Intended users are health care professionals only and are defined by the application of the compatible stereotactic system.

• Intended patients

Intended patients are defined by the compatible stereotactic system according to its application and are not given by the needle guide bushing itself.

Description

The needle guide bushing STBG-CY is made of medical-grade stainless steel.
The needle guide bushings are sterile, packed in pairs in peel pouches and are available in various diameters (see Table 1).
Sales units are additionally packaged in protective cardboard box.

Compatibility

The needle guide bushing STBG-CY is specially designed for use with the compatible Cytoguide stereotactic system.
For compatible needles and cannulas, see Tab. 1.

Application

1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. Further instructions should be followed according to the Cytoguide stereotactic system instructions for use.

2. Preparation of the medical device

Prior to use, the gauge size (outer diameter) required for the application to be performed must be determined and the compatibility of the needle guide bushing with the needle or cannula used must be ensured. See Tab. 1.

3. Removal from sterile package

Check sterile packaging and the product for expiration date and visual damages before use.
Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly  and remove the product aseptically.

4. Application of medical device

Insert the two needle guide bushings into the receptacles provided at the Cytoguide stereotactic system.
Insert the compatible needle or cannula through the opening of the needle guide bushings and apply according to instructions for use.
All application steps must be carried out in accordance with the instructions for use of the Cytoguide stereotactic system and the compatible medical devices.

5. Removal and disposal

After procedure remove needle guide bushings from stereotactic system and dispose of in accordance with local provisions and safety regulations.

REACH declaration at <https://www.bipmedical.com/en/about-us/regulations/> must be observed.

Warnings

- The needle guide bushing STBG CY may only be used in combination with the Cytoguide stereotactic system. A combination with other stereotactic systems is not permitted and can lead to serious injury to the patient.
- The needle guide bushing must not be used and must be disposed of if the product is damaged or if the expiration date has elapsed or if the sterile packaging is damaged or unintentionally opened before use.
- The needle guide bushing is intended for single use only and must not be reused. Reuse and re-sterilization impair the safe functionality and can lead to cross-infection of patient or material failure of the product.
- The needle guide bushing is not MR-compatible.
- The used needle guide bushing must be disposed of with suitable measures to prevent injuries or biohazardous contaminations.

Needle Guide Bushing STBG-CY

Precautions

- The warranty applies exclusively to the proper intended use of the product.
- Please also note the operating instructions of compatible products. The manufacturer assumes no liability and warranty for use with non-compatible medical devices.
- Before use inspect each product for damages that would prevent proper functionality. Do not use damaged products.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

Environmental storage conditions

Temperature		10 – 25°C
Humidity		20 – 50% r.h (non-condensing)
Shelf-life		3 years

NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

Technical data / available variants

See Tab. 1.

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©2025 BIP GmbH, Germany, all rights reserved.